

Shoprider TE888IX


enjoy the freedom...



This is our newest addition to the Shoprider range of scooters. A midsize scooter full of features and clothed in a new modern body with integrated replaceable bumpers and silver spoked mag wheels fitted with puncture-less tyres. A powerful new LED headlight, new indicators and tail lights as well as a new engine and tiller handlebar design complete the package.

FEATURES

- Powerful 1.0hp motor
- Headlights & reflectors for night use
- Modern ergonomically designed steering tiller
- On/Off security key
- Large range of accessories
- 4 wheel suspension
- Available in Red or Blue



RED BLUE



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



SPECIFICATIONS

width	600mm
length	1250mm
height (headrest lowered)	1200mm
weight (incl battery)	87kg
speed*	0-10 kph
range*	30 km
max. climb rate*	9°
turning radius	1420mm
max. load	136 kg

All specifications are correct as at March 2016 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.