

Shoprider TE888NR

enjoy the freedom...



Stability, simplicity of operation and affordability are just a few of the features on one of our most popular models.

This midsize scooter with its fully adjustable seat and easy to operate fingertip controls has a range of up to 30 kilometers making short or long journeys a smooth comfortable ride.



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



Shoprider New Zealand Ltd • P.O. Box 35-059 Shirley • Christchurch 9 • New Zealand

FEATURES

- Speed dial allows you to alter top speed
- On / Off security key
- Battery level indicator with charger socket mounted high on tiller for easy access
- Horn
- Shopping basket
- New design adjustable steering tiller
- Powerful headlight for added security
- Deluxe adjustable seat - rotates, slides backward / forward with folding back and lift up armrests

SPECIFICATIONS

width	550 mm
length	1250 mm
weight (incl battery)	82 kg
speed*	0-8 kph
range*	30 km
max. climb rate*	9°
max. load	136 kg

All specifications are correct as at March 2016 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.