

Shoprider TE889SL

enjoy the freedom...



This wonderfully comfortable scooter owes its super smooth ride to its 4 wheel suspension its larger 330mm tyres and the adjustable captains chair.

It also comes equipped with Delta Tiller handlebars for comfortable hands while the heavy duty, powerful 1.0hp motor was made for very demanding use.

FEATURES

- Adjustable captains chair with head rest and fold up arm rests
- Headlights & reflectors for night use
- Modern ergonomically designed steering tiller
- On/Off security key
- Large range of accessories
- Independent suspension with Shock Absorbers

SPECIFICATIONS

width	640 mm
length	1300 mm
height	1300 mm
weight (incl battery)	108 kg
speed*	0-12 kph
range*	38 km
max. climb rate*	9°
turning radius	1500 mm
max. load	136 kg

All specifications are correct as at March 2016 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER®
New Zealand

Authorised Dealer

