

Shoprider TE889XLSN

enjoy the freedom...



Built for demanding use the TE889XLSN can carry up to 225kg.

The leatherette pillow-top captains seat provides a luxurious ride while the 4 wheel suspension and large wheels offer greater ground clearance and a degree of comfort that is difficult to match.

FEATURES

- Deluxe adjustable reclining seat
- Headlights, tail lights, hazard lights & indicators
- Modern ergonomical steering tiller
- Horn
- 4 wheel suspension
- Hi / Lo ratio speeds
- Separate charger
- Strong 2hp motor

SPECIFICATIONS

width	690 mm
length	1470 mm
height	1300 mm
weight (incl battery)	147 kg
speed*	0 - 15 kph
range*	54 km (75amp batteries)
safe climb angle*	12°
turning radius	1829 mm
max. load	225 kg

All specifications are correct as at March 2016 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



Shoprider New Zealand Ltd • P.O. Box 35-059 Shirley • Christchurch 9 • New Zealand