

Shoprider TE9

enjoy the freedom...



This is our newest addition to the Shoprider range of scooters.

A midsize scooter full of features such as a high / low speed switch as well as a speed dial allowing the top speed to be altered between 0 - 12.5 kph. This new scooter has the same powerful headlight as the TE888NR combined with a new high back seat. It also comes with throttle lever extensions.



Sign up to Shoprider's AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



FEATURES

- Gas tiller for easy angle adjustment **NEW**
- 5 stage adjustable rear suspension **NEW**
- Front and rear suspension
- Two rear view mirrors for extra visibility
- High back captain seat with adjustable headrest and armrests
- Rounded handlebar to make driving easy and comfortable
- Front headlight, tail lights and indicators

SPECIFICATIONS

width	700 mm
length	1270 mm
height	1270 mm
total weight	103 kg
speed*	0 - 12.5 kph
range*	40 km
max. climb rate*	9°
turning radius	1500 mm
max. load	159 kg

All specifications are correct as at March 2016 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.