

# Shoprider TEGK9

enjoy the freedom...



This stylish 4 wheel scooter gets around those tight indoor areas easily. It manoeuvres easily inside your house getting through the narrowest of doors.

This new model from Shoprider is very quiet.

Its compact lightweight design also allows it to be dismantled easily for transporting.



Sign up to Shoprider's AA Mobility Care programme and make sure you Get Home Safely with the AA

**SHOPRIDER**<sup>®</sup>  
New Zealand

Authorised Dealer



## FEATURES

- Speed dial allows you to alter top speed
- On/off security key
- Battery level gauge
- Non puncture tyres
- Adjustable armrests
- Easily disassembled
- Choice of colours

## SPECIFICATIONS

width	560 mm
length	1050 mm
height	880 mm
total weight	56 kg
speed*	0-6 kph
range*	20 km
max. climb rate*	9°
turning radius	1150 mm
max. load	136 kg

All specifications are correct as at March 2016 but are subject to change without notice. \*Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.