

Shoprider TE-779XLS

enjoy the freedom...



This Deluxe 3 wheel provides the comfort of a large scooter with the added advantage of the deluxe suspension system.


The 3 wheel system is very maneuverable and easy to drive and has a weight capacity of 185 kg combined with a high back with headrest sliding/rotating/height adjustable seat for user comfort.



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



FEATURES

- Powerful 1.0hp motor
- LED head-light and rear lights, indicators and hazard lights
- High, medium and low speed settings allows you to alter top speed
- On/Off security key
- Large range of accessories
- Front and rear suspension

SPECIFICATIONS

width	680mm
length	1525mm
height (headrest lowered)	1300mm
weight (incl battery)	140kg
speed*	0-15kph
range*	65km
max. climb rate*	9°
turning radius	1650mm
max. load	185kg

All specifications are correct as at December 2022 but are subject to change without notice. *Maximum speed, range and climb are subject to terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown.